

How it works


MRC has pioneered the use of appliances to correct myofunctional habits in growing children and has proven success in orthodontic correction without braces. This treatment can also lead to better facial development in growing children. The key to this treatment is correcting the position and function of the tongue, obtaining correct nose breathing and retraining the oral muscles to function correctly.

Myobrace® appliances effectively train the tongue to position correctly in the upper jaw, retrain oral musculature and exert light forces to align the teeth.

Goals of treatment:

- Lips together at all times except when speaking or eating.
- Breathe through the nose, to assist the development of the upper and lower jaws, and to achieve the correct bite.
- No lip activity when swallowing, which allows the front teeth to develop correctly.
- Improved dental alignment.
- Improved facial development.

For the Myobrace® to work, you must correct your poor myofunctional habits.



Other Myobrace® systems:



There are Myobrace® systems suitable for a range of patient ages and needs.

- **Myobrace® for Juniors** ages 3 - 5
- **Myobrace® for Teens** ages 10 - 15
- **Myobrace® for Adults** ages 15+

Talk to your doctor today about Myobrace® or visit www.myobrace.com.

myobrace®
SUPREME
DENTAL
CONCEPT

FOR KIDS

**Straight teeth
the natural way.**

Early treatment
can prevent the
need for braces.

AGES
5-10

myobrace®

For more information, visit www.myobrace.com
MYOBACE® BY MYOFUNCTIONAL RESEARCH CO.


www.myobrace.com

What causes orthodontic problems?

From 6-10 years of age, children are getting their first permanent teeth. At this stage the teeth can become crowded and the jaws misaligned.

Habits such as mouth breathing, incorrect lip and tongue function (known as poor myofunctional habits) are the major causes of poor jaw development and tooth alignment. Research shows that these problems continue through to the permanent dentition, requiring complex orthodontic treatment. Braces and extractions do not resolve these problems and are best treated as soon as they are evident.

These habits are easier to correct in the growing years. Once the causes are corrected natural growth and dental alignment can take place.



Is orthodontics with braces the only treatment option?

Since 1989 Myofunctional Research Co. (MRC) have developed appliances that provide a more effective option for the growing child with orthodontic problems. The appliances are designed not only to straighten teeth, but also to treat the underlying causes of crooked teeth by correcting poor oral habits such as mouth breathing, tongue thrust and incorrect swallowing patterns.

Unlike traditional orthodontics, the primary goal of MRC's approach to treatment is to achieve naturally straight teeth for life – often without the need for braces or retainers. Treating the causes of crooked teeth also improves general health and allows children to reach their maximum genetic potential.

Myobrace® for Kids

Myobrace® for Kids is a three-stage appliance system designed specifically to correct poor oral habits while treating upper and lower jaw development problems. *Myobrace® for Kids* is most effective after a child's permanent front teeth have come through and before all the permanent teeth have erupted.

The appliances correct poor oral habits which are the cause of crooked teeth, while providing light forces to help expand the upper jaw and align the teeth into their correct natural position.

Treatment considerations

Treatment time can differ depending on the biological adaptability of each patient. For a successful outcome the *Myobrace®* appliance must be used daily for 1-2 hours plus overnight while sleeping.

Regular monitoring by your dentist or orthodontist is essential for successful treatment and to avoid unwanted treatment complications.

Normal vs Crowded teeth

The function of the tongue, cheeks and lips determines tooth position. Correct facial growth depends on correct function and breathing patterns.



CROWDED
UPPER ARCH

NORMAL
UPPER ARCH

NOTE: TOOTH SIZE HAS NOT CHANGED

Myobrace for Kids design features



This patient's spacing was treated successfully using *Myobrace®* treatment methods.